

RetrieveRx®

Eliminate the clutter and guesswork of storing and locating filled prescriptions.

Integrated with QS/1's NRx® and PrimeCare® Pharmacy Management Systems, RetrieveRx electronically organizes prescriptions and assigns them to bags for quick, easy retrieval. As with other QS/1 products, you have the flexibility to modify RetrieveRx to fit your store's needs.

Integrated System Saves Time

As an integrated part of QS/1's Pharmacy Management Systems, RetrieveRx is easy to learn and use. It also simplifies service and support since there is only one call to make. Integration also eliminates time-consuming redundancies that occur between systems from different vendors. It also aids your customer service by enabling your staff to answer price and pick-up questions right from their workstations.

One Click Return to Stock

Integration also means a simplified return to stock. The system lights the prescription bags that need to be returned to stock. One mouse click in the pharmacy management system voids a transaction and reverses actions related to the transaction including claims, workflow records, IVR callouts and inventory.

Sturdy Hardware that Lasts

The hardware utilized is the RxMedic® ARS® (Automatic Retrieval System), a wireless retrieval system. The hardware includes heavy-duty bags and handles. The crystal clear bags remain closed while on the rack and light up whether on the hanger or not.



- Helps pharmacies store, then locate, filled prescriptions.
- Integrated with QS/1's NRx and PrimeCare Pharmacy Management Systems.
- Handles individual or batched prescriptions.
- Accounts for refrigerated prescriptions.
- One click return to stock.
- Eliminates redundancies.
- Reduces patient wait.
- Enables staff to answer price and pick-up questions from their workstations.
- Durable wireless hardware system.

Rx Summary: #06002163 for TEST, PATIENT

Verify Prescription

NDC: 00186-5020-31
Scan NDC: 00186-5020-31 [Complete] [Close]

Rx Number: 06002163

Patient: TEST, PATIENT (old Folks Home)
Charge: TEST, PATIENT
Prescriber: TEST, DOCTOR
Primary: TEST PLAN
USC Plan: USC
Secondary:
Tertiary:

Drug: NEXIUM 20 MG CAPSULE Expiration: 02/20/2013
Generic: []
Sig: TAKE ONE CAPSULE DAILY.
HOA: [] Frequency: [] Cycle Fill: []
Rx Notes:

Assign Bag

Scan Bag: [] [Find] [Continue]

Other Information

Auth Number: [] Status: Active
DAW: Y 0 - No Product Special: []

Pricing Information

Price: 219.17 [Full] [Part]
Acquisition: \$112.45 Secondary: \$0.00
Profit %: 94.91% Tertiary: \$0.00
Copy: \$0.00 Third Party: \$219.17

Drug Description

This drug is a amethyst, oblong-shaped, capsule imprinted with NEXIUM 20mg on one side.

Additional Rx Information

Place of Service: []
Type of Service: []
Special Pkg Ind: 2 - Manufactur
Rx Origin: 1 - Written

Messages

No Refills Remaining
SP Field Is Blank
Prescriber Phone Number - 608-647-6161
Warning - Early Refill 30 Days
Charge Limit Exceeded

document management interface error - document system is not running.

Advantages of Integrated System

- Easier to learn and use because it is integrated with the pharmacy management system.
- Allows staff to answer customer questions from their workstations.
- Eliminates redundancies between will call and pharmacy management systems.
- Makes return to stock simple with one click of the mouse.
- Low battery status.
- Include multiple prescriptions in bags.
- Include prescriptions for multiple family members.
- Bags are detected automatically when scanned into the pharmacy system, requiring no prior setup for the system to recognize the bag.
- 35 different lighting configurations that can identify which bag is for which operator.

Hardware Specifications

- Handles and bags are made of heavy-duty plastic.
- Two bag sizes: 12"x 9" and 12"x16"
- Wireless infrared transmitter.
- Battery life is 18-24 months.
- Easy-open handle that stays closed while on the rack.
- Uses disposable 3V coin batteries.
- Can assign 10 bags per 1 foot section.
- Hanging rods available in 2', 3' and 4' lengths.



800.231.7776 www.qs1.com